
295 Cartridge Box
UniDose & Quickfill
ref 1110

EN: Instructions for use
NL: Gebruiksaanwijzing
IT: Istruzioni per l'uso

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Date : 01/2025
 DR/0490201 V6

CE 2797

INTENDED USE
The Cartridge (UDI 05430002956008) is a single-use class I sterile medical device intended to contain radiopharmaceutical drugs and to administer its contents to humans through an injection line.

INTENDED POPULATION
Medicine nuclear patients as identified by healthcare professionals.

INTENDED USERS
The Cartridge is intended to be used by qualified pharmaceutical and medical personnel specifically trained in its use and in handling radioactive drugs.

CONTRAINDICATIONS
None known for the devices concerned.

ADVERSE EVENTS/COMPLICATIONS
None known for the devices concerned.

WARNINGS/CAUTION
Do not use the product if its packaging is damaged and/or opened. Single-use device; do not reuse or resterilise. This could compromise the integrity of the device, cause malfunctions and/or infections or injuries.
When injecting radioactive substances, precautions must be taken to protect the personnel and the public against radiation (shielding). Use with the compatible injection line as specified.

INSTRUCTIONS FOR USE
1. Insert the 295 cartridges box into the dispenser as described in the dispenser manual.
2. After collection of a single radiopharmaceutical dose in a cartridge connect to the injection line according to the dispenser instructions and proceed with injection into the patient.
3. Discard in approved collector after use.
Note: the letters ID on the box label identify the box number in its batch; V3 is the software version of the RFID-TAG (see Unidose manual).
Users should report any serious incident related to the device to the Manufacturer and National Competent Authority.

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